

Sample Grinder DF-01A



DF-01A is our self-developed supporting device for spectrometers. It's designed based on sample quantity and material and lab application. Several models are available now.

- 1) Grinding plate: 350mm in diameter attached by single or double trays, sandpaper or sanding disc.
- 2) Single or double grinding plates. You can attach abrasive disc or sandpaper to the double plates by need.
- 3) Motor speed: 1400r.p.m, electromagnetic brake
- 4) Max power: single tray 2.7 kw; double trays 4.2 kw.
- 5) Self-cleaning function
- 6) Hydraulic hood lifting
- 7) Motor will be off automatically when the hood is lifted. It will also be protected from over electrical current caused by short circuit.
- 8) 900*530*880(mm), 100kg.